A one-touch solution for deep vein thrombosis



Deep vein thrombosis (DVT) occurs when a blood clot forms in a deep vein, often in the legs. Some of it may travel through the veins to the lungs and block its blood supply, resulting in a pulmonary embolism or pulmonary hypertension, both of which can be fatal.

By applying intermittent pressure through inflatable cuffs, the Landing IPC System creates a pressure gradient in the limb. This prevents stasis, accelerates venous blood flow, and enhances fibrinolysis. Other features include:

- All-in-one controller with intuitive display
- One-touch operation
- Adjustable pressure
- Automated alarm system
- Quiet operation





ELANDING QOMO

Landing Qomo 46950 Magellan Drive Wixom, MI 48393

USA

Phone: +1-866-990-7666 ext. 3147

Fax: +1-248-281-1998

Email: info@landingqomo.com Web: www.landingqomo.com

